



# Single Use Circular Stapler

## Patented Cutting Blade, Breakthrough Innovation in Double Stapling Technology

- Unique staple cutting technology, no sticking of staples after firing
- Safer in low rectal sphincter-preserving procedure
- No damage to anastomosis



### Major Features

- Applicable to end-to-end, end-to-side anastomosis of GI tract and esophagogastric devascularization procedure
- Unique 23mm size, most appropriate for esophageal supra-arch anastomosis
- Precise adjusting of closed staple height, applicable to all tissues
- Gastric tube port in 23 and 26mm models to simply gastric procedures

**Golden Stapler Surgical Co., Ltd.**



*Advancing Surgical Technology  
to Benefit More People*

### Patented Staple Cutting Technology

- No sticking of staples after firing
- No damage to anastomosis

### 5.3mm Open Leg Length

- Staple formation height range: 1.0-2.5mm
- Satisfy stapling of various tissue thickness of gastroenterostomy

### Thin and Flat anvil head

- Easy to insert into enteric cavity & remove from anastomosis
- Avoid damage to intestinal wall & anastomosis

### Short Distance between Firing and Grabbing Handles

- Simple firing mechanism no exertional effort required

### Silicon Coated Handle

- Ergonomic design of great hand touch feeling



### Gastric Tube Port

- Load gastric tube in advance, pull it out together with the circular stapler after surgery
- Simplified procedure, reduced duration

### 5 Different Sizes

- Establish the most appropriate anastomosis
- Allow for various gastroenterostomy

### Conjoined Trocar

- Simplify surgical procedures
- Facilitate the lower anterior resection

### Plastic Cutting Ring






- Tactile & audible responses
- Indication of safe firing sequence completion

### Magnified Green Indication Window with a Gap Setting Scale

- Controllable tissue compression and staple formation height
- Effectively achieve suturing & homeostasis purposes

### Reverse Conjoined Trocar

- Applicable for double stapling technique

Code	Size (MM)	Open Staple Height (MM)	Staple Dia. (MM)	Adjustable Closed Staple Height (MM)	Colour	Staple Qty	Circular Knife Dia. (MM)	Staple Rows
CSB21	21	5.3	0.28	1.0~2.5	Jacinth 	18	13	2
CSB23	23	5.3	0.28	1.0~2.5	Orange 	18	15	2
CSB26	26	5.3	0.28	1.0~2.5	White 	20	17	2
CSB29	29	5.3	0.28	1.0~2.5	Blue 	24	20	2
CSB32	32	5.3	0.28	1.0~2.5	Green 	30	22	2

\* Staple from titanium material

**Golden Stapler Surgical Co., Ltd.**

Add: Building 7A, Jiangsu Wujin Sci-Tech Pioneer Park, 256 Mid Mingxin Road, Hutang Town, Wujin District, Changzhou, Jiangsu Province, China 213164  
 Tel: +86 519 85070251 Fax: +86 519 85070252  
 Web: <http://www.goldenstapler.com>