




# Single Use Hemorrhoidal Circular Stapler

**Multi-functions Design, Tissue Stopper  
32mm Lower than PPH03(4mm)**

- Designed for both hemorrhoids and STARR procedures
- Deepest and largest staple housing, easier to contain more cut tissue



## Major Features

- Accommodate to II~IV internal hemorrhoids
- Accommodate to STARR (Stapled Transanal Rectal Resection) procedure to treat Obstructed Defecation Syndrome (ODS)
- Protect anal cushion and retain precise faeces control
- Simple and minimally invasive to use
- Faster procedure, less pain, quicker recovery

**Golden Stapler Surgical Co., Ltd.**



*Advancing Surgical Technology  
to Benefit More People*

# Single Use Hemorrhoidal Circular Stapler

## Staple Housing is 6mm Longer than PPH03

- Staple housing size is larger than PPH03 over 30%
- Staple housing is able to contain much more cut tissue

## Plastic Cutting Ring

- Optimized tactical & audible responses during firing

## Transparent Anoscope & Butterfly Anal Dilator

- Visible operation and position selection during purse-string procedure
- Butterfly shaped dilator, applicable to obese patient or strong hip muscle

## Butterfly Shaped Nut

- Decent looking, better touch feeling, easier to operate
- Precise adjusting of closed staple height

## Magnified Green Indication Window with a Gap Setting Scale

- Controllable staple formation height between 0.8-2.0 mm
- Accommodate to PPH and STARR procedures
- Effectively achieve suturing and homeostasis purposes

## Tissue Stopper Design

- 32mm Lower than PPH03 (4mm)
- Much easier to retain more cutting tissue into staple housing

## Less Distance between Firing and Grabbing Handle

- Facilitate female surgeons' easier operation
- Silicon handle avoids slippage and maintains a great hand touch feeling

## Larger Suture Thread Ports

- Significantly facilitate dragging cutting tissue into staple housing

## Smooth Design of Staple Housing Edge & More Staples

- Drag cut tissue into staple housing without damage
- 32mm model contains 4 more staples than PPH03, better hemostasis

## Non-detachable Central Rod Design

- One-piece central rod design to avoid anvil fall off
- Solid and non-detachable anvil offers premium staple formation basis



Code/Model	Open Staple Height (MM)	Staple Qty	Closed Staple Height (MM)	Outer Dia. (MM)	Circular Knife Dia. (MM)	Staple Rows	Remark: PPHB and PPHH models also available
PPHD32	4.5	32	0.8-2.0	32.5	24	2	PPHD: Butterfly Anal Dilator; Opaque Staple Housing
PPHD34	4.5	34	0.8-2.0	34.5	26	2	PPHB: Round Anal Dilator; Opaque Staple Housing
PPHD36	4.5	36	0.8-2.0	36.5	28	2	PPHH: Butterfly Anal Dilator; Transparent Staple Housing

\* Staple from titanium material

**Golden Stapler Surgical Co., Ltd.**

Add: Building 7A, Jiangsu Wujin Sci-Tech Pioneer Park, 256 Mid Mingxin Road, Hutang Town, Wujin District, Changzhou, Jiangsu Province, China 213164  
 Tel: +86 519 85070251 Fax: +86 519 85070252  
 Web: <http://www.goldenstapler.com>