

Single Patient Use Endo Cutter and Reload

Innovative Stapling Technology: **Stepped Closed Staple Heights with Unequal Height Differences**

Good hemostasis in the tissue cutting edge; better blood-supply in the healthy tissue

Major Features

- Applicable to the cutting and stapling of GI tract, lung, vessel or mesentery in open and endoscopic surgeries
- Available of 60, 160 and 260mm hand pieces; 45 and 60mm reloads of all 5 colors;
 30mm white reload with curved tip. Full range of choice
- o Any hand piece and reload are mutually compatible. Reduce storage cost and simplify loading of reload
- o Reload tip to articulation joint is 2mm shorter than traditional products. Easier access to deeper pelvic cavity
- o Fixed anvil instead of fixed cartridge. Easier placement of of device in target position

Golden Stapler Surgical Co., Ltd.

Advancing Surgical Technology to Benefit More People











Single Patient Use Endo Cutter and Reload

Cutting Line Output Output

Unequally Stepped Closed Staple Heights

- · Increasing closed staple heights from cutting edge to outer side
- · Unequal height differences, best balancing between hemostasis and blood-supply

Curved Tip Design

- · Unique curved tip in 30mm white reload
- · Much easier to position reload into narrow tissue space



8.5cm in GoldenStapler 45mm Reload

"Wrist"Instead of "Elbow" & Fixed Anvil

- GoldenStapler Reload's joint in "wrist" position instead of traditional "elbow" position
- Distance between reload tip and joint 2cm shorter, capable of entering the deepest pelvic position
- Fixed anvil design instead of fixed cartridge, easier to place under target tissue



Full Portfolio of Reloads

- Grey, white, blue, yellow and green for both 45 and 60mm reloads; Unique 30mm white reload with curved tio
- · All reloads compatible with any hand piece
- Articulated up to 45°left and right, easier to be placed into target position
- · Non-articulated reload also available
- Integrated cutting knife, no contamination during firing. Comply with tumour-free principle

Fine Designed Reload

- Reload can clamp the tissue without damage
- Green reload can go through 12mm endoscopic port, less traumatic
- Safety mechanism for empty reload to prevent misfiring
- Stapling line longer than cutting line, zero bleeding after firing

Universal Hand Pieces of Premium Quality

- 60, 160 and 260mm models for different situations
- 60mm hand piece can be applied in open surgeries;
 160mm hand piece in common endo-surgeries;
 260mm hand piece for obese patients
- · All models compatible with any reloads
- Shaft can be rotated 360°, easy to find the most convenient angle for firing

Hand Piece:

Code	Shaft Length (MM)	Remark	
EC-60-I	60		
EC-160-I	160		
EC-260-I	260		
EC-60	60		
EC-160	160	Customer Made	
EC-260	260		

Reload:

Code	Articulated Reload Colour	Staple Rows	Open Staple Height (mm)	Closed Staple Height (mm)	Remark
ECRA45M	Grey	6	2.0	0.8 - 0.8 - 0.8	Customer Made
ECRA60M		6			Custoffier Made
ECRC30T	White	4	2.5	0.8-0.8	Curved Tip
ECRC30W		6		0.8 - 1.0 - 1.1	Curved Tip
ECRA45W		6			
ECRC45W		6			Curved Tip
ECRA60W		6			
ECRC30B	Blue	6	3.5	1.3 - 1.5 - 1.6	Curved Tip
ECRA45B		6			
ECRC45B		6			Curved Tip
ECRA60B		6			
ECRC60B		6			Curved Tip
ECRA45Y	Yellow 🦲	6	3.8	1.6 - 1.8 - 1.9	
ECRA60Y		6			
ECRA45G	Green	6	4.1	1.8 - 2.0 - 2.1	
ECRA60G		6			

^{*} Non-articulated Reloads ECRB codes also available

^{*} Staple from titanium material